


Treff Vial HP

(Transport and Autosampler Vial)





CleanRoom®
PURE

The Treff Vial HP with septum cap is specially designed to meet required parameters of modern transport and autosampler vials.

General Characteristics – Treff Vial HP

Thanks to its special design and 2K material combination it is capable of withstanding up to 5 bar internal pressure and has a 2 bar resealing tightness after needle penetration of the septum. Its outstanding DMSO tightness makes the Treff Vial HP to an ideal storage and transport vessel for a wide range of research, development, quality control and production laboratories. Each vial features an individual unique 2D tube number to ensure retraceability at any time.

All TreffLab Vial HP feature the following characteristics:

- › Outstanding dynamic transport tightness and resealing properties
- › Autosampler and transport system approved
- › 100% virgin medical-grade PP only, no additives, providing high chemical quality
- › Microbiological controlled production in cleanrooms class 8, ISO 14644-1
- › Free of RNase, DNase, human DNA, endotoxins, PCR inhibitors and ATP
- › Capacity: up to 350 µl
- › Up to 5 bar / 72.5 psi internal pressure tightness
- › Resealing tightness after needle penetration of septum: 2 bar / 29 psi
- › Temperature range: - 80°C up to 100°C
- › Ready to Use (CleanRoom Pure®)

Treff Vial HP 350 µl *transport and autosampler vial*

- ▶ Leak safe sample transport vessel for fully automated delivery systems
- ▶ Individual unique 2D barcode



1

Product Description	Material	Color	Quantity	Order No.	Image No.
Treff Vial HP 350 µl	PP	clear	5'000 10 bags of 500 Tubes 10 bags of 500 Caps	96.11418.9.00	1



Nolato Treff AG

Taastrasse 16, 9113 Degersheim

Switzerland

Phone: +41 71 372 55 55

E-Mail: treff.ag@nolato.com, trefflab@nolato.com

www.nolato.com/treff

www.trefflab.com

