


Technical Drawing

96-Well PCR Plate 0.15 ml, PP,
low profile, full skirt, clear

Order number 96.11203.9.01




L	127,8 mm	Yn	74,3 mm
W	85,5 mm	P	9,0 mm
H	15,5 mm	D	5,5 mm
X1	14,4 mm		
Xn	113,4 mm	C	12
Y1	11,2 mm	R	8

Dimensions : Millimeters

Technical Drawing

96-Well qPCR Plate 0.15 ml, PP, low profile, full skirt, white

Order number 96.11197.9.01




L	127,8 mm	Yn	74,3 mm
W	85,5 mm	P	9,0 mm
H	15,5 mm	D	5,5 mm
X1	14,4 mm		
Xn	113,4 mm	C	12
Y1	11,2 mm	R	8

Dimensions : Millimeters

Technical Drawing

96-Well PCR Plate 0.2 ml, PP,
 standard profile, half skirt, clear
 Order number 96.11204.9.01




L	127,8 mm	Yn	74,2 mm
W	85,5 mm	P	9,0 mm
H	20,7 mm	D	5,2 mm
X1	14,4 mm		
Xn	113,4 mm	C	12
Y1	11,2 mm	R	8

Dimensions : Millimeters

Technical Drawing

96-Well qPCR Plate 0.2 ml, PP,
standard profile, half skirt, white

Order number 96.11198.9.01



L	127,8 mm	Yn	74,2 mm
W	85,5 mm	P	9,0 mm
H	20,7 mm	D	5,2 mm
X1	14,4 mm		
Xn	113,4 mm	C	12
Y1	11,2 mm	R	8


Dimensions : Millimeters

Technical Drawing

96-Well PCR Plate 0.15 ml, PP,

low profile, half skirt, clear

Order number 96.11205.9.01




L	127,8 mm	Yn	74,2 mm
W	85,5 mm	P	9,0 mm
H	15,3 mm	D	5,2 mm
X1	14,4 mm		
Xn	114,4 mm	C	12
Y1	11,2 mm	R	8

Dimensions : Millimeters

Technical Drawing

96-Well qPCR Plate 0.15 ml, PP, low profile, half skirt, white

Order number 96.11199.9.01




L	127,8 mm	Yn	74,2 mm
W	85,5 mm	P	9,0 mm
H	15,3 mm	D	5,2 mm
X1	14,4 mm		
Xn	114,4 mm	C	12
Y1	11,2 mm	R	8

Dimensions : Millimeters

Technical Drawing

96-Well PCR Plate 0.15 ml, PP,
 low profile, half skirt, elevated rim, clear
 Order number 96.11206.9.01




L	127,2 mm	Yn	74,1 mm
W	85,2 mm	P	9,0 mm
H	15,3 mm	D	5,2 mm
X1	14,1 mm		
Xn	113,1 mm	C	12
Y1	11,1 mm	R	8

Dimensions : Millimeters

Technical Drawing

96-Well qPCR Plate 0.15 ml, PP,
 low profile, half skirt, elevated rim, white
 Order number 96.11200.9.01




L	127,2 mm	Yn	74,1 mm
W	85,2 mm	P	9,0 mm
H	15,3 mm	D	5,2 mm
X1	14,1 mm		
Xn	113,1 mm	C	12
Y1	11,1 mm	R	8

Dimensions : Millimeters

Technical Drawing

96-Well PCR Plate 0.2 ml, PP,
standard profile, no skirt, clear

Order number 96.11207.9.01




L	120 mm	Yn	71,5 mm
W	80 mm	P	9,0 mm
H	20,7 mm	D	5,2 mm
X1	97,5 mm		
Xn	109,5 mm	C	12
Y1	63,5 mm	R	8

Dimensions : Millimeters

Technical Drawing

96-Well qPCR Plate 0.2 ml, PP,
standard profile, no skirt, white

Order number 96.11201.9.01




L	120 mm	Yn	71,5 mm
W	80 mm	P	9,0 mm
H	20,7 mm	D	5,2 mm
X1	97,5 mm		
Xn	109,5 mm	C	12
Y1	63,5 mm	R	8

Dimensions : Millimeters

Technical Drawing

96-Well PCR Plate 0.15 ml, PP, low profile, no skirt, clear

Order number 96.11208.9.01





L	120 mm	Yn	71,5 mm
W	80 mm	P	9,00 mm
H	15,3 mm	D	5,2 mm
X1	10,5 mm		
Xn	109,5 mm	C	12
Y1	8,5 mm	R	8

Dimensions : Millimeters

Technical Drawing

96-Well qPCR Plate 0.15 ml, PP, low profile, no skirt, white

Order number 96.11202.9.01



L	120 mm	Yn	71,5 mm
W	80 mm	P	9,00 mm
H	15,3 mm	D	5,2 mm
X1	10,5 mm		
Xn	109,5 mm	C	12
Y1	8,5 mm	R	8



Dimensions : Millimeters

Technical Drawing

96-Well qPCR Plate 0.15 ml, PP,


low profile, half skirt, white (for Roche® LightCycler® 480)

Order number 96.10867.9.01, order number 96.10868.9.01



Dimensions : Millimeters

Technical Drawing
384-Well PCR Plate 0.40 µl, PP,
full skirt, clear
 Order number 96.09974.9.01




L	127,76 mm	Yn	76,49 mm
W	85,48 mm	P	4,50 mm
H	9,50 mm	D	3,20 mm
X1	12,13 mm		
Xn	115,63 mm	C	24
Y1	8,99 mm	R	16

Dimensions : Millimeters

Technical Drawing

384-Well qPCR Plate 0.40 µl, PP,
 full skirt, white (for Roche® LightCycler® 480)
 Order number 96.10180.9.01



L	127,76 mm	Yn	76,49 mm
W	85,48 mm	P	4,50 mm
H	9,50 mm	D	3,20 mm
X1	12,13 mm		
Xn	115,63 mm	C	24
Y1	8,99 mm	R	16

Dimensions : Millimeters